



厦门为正生物科技股份有限公司
XIAMEN WIZ BIOTECH CO., LTD

Declaration on the detection validity of omicron variant

To Who It May Concern,

For the newly discovered variant of novel coronavirus Omicron (B.1.1.529), we have analyzed WIZ SARS-CoV-2 Antigen Rapid Test Kit and came to a conclusion: there is no off-target or missing detection occurred in the latest test for variant of Omicron (B.1.1.529). WIZ SARS-CoV-2 antigen rapid test adopts double antibody sandwich method to detect Nucleocapsid protein of SARS-CoV-2 viruses, Omicron NP protein mutation sites were P13L, 31-33, R203K, G204R. WIZ SARS-CoV-2 antigen rapid test detectable sites are not at these mutation sites. Our team will continue to track the latest variant of novel coronavirus to ensure stable supply.

XIAMEN WIZ BIOTECH CO., LTD.

2021.11.30

